

MOLECULAR DETECTION OF *BACILLUS ALTITUDINIS* 19RS3 AND T5S-T4 ISOLATED FROM YERBA MATE (*ILEX PARAGUARIENSIS* ST. HIL.) USING STRAIN-SPECIFIC PRIMERS

I Simpósio de Microbiologia de Rondônia: Saúde, Ambiente e Inovação., 1ª edição, de 23/03/2021 a 25/03/2021
ISBN dos Anais: 978-65-86861-91-4

CORTESE; ILIANA JULIETA ¹, CASTRILLO; MARIA LORENA ², ONETTO; ANDREA LILIANA ³, ZAPATA; PEDRO DARIO ⁴, LACZESKI; MARGARITA ESTER ⁵

RESUMO

Introduction: The genus *Bacillus* presents a great diversity of species that are widely distributed in the environment. It is one of the most studied genera and was shown to improve plant growth through a combination of mechanisms. It is used as a biofertilizer that holds promise to make sustainable agricultural practices and ecologically safe. In our previous studies, two endophytic endospore-forming bacteria coded as *Bacillus altitudinis* 19RS3 and T5S-T4 were isolated from *Ilex paraguariensis* St. Hil roots and selected for their plant growth-promoting (PGP) properties *in vitro* and *in vivo*. Strain-specific primers to detect *B. altitudinis* 19RS3 and T5S-T4 strains were designed by whole-genome analysis. **Objective:** This study aimed to test the pre-designed strain-specific primers with DNA isolated from *B. altitudinis* 19RS3 and T5S-T4. **Methods:** Genomic DNA from liquid cultures incubated at 30°C for 24 h was extracted by using Sambrook work protocol modified. For the molecular detection, the primers designed for 19RS3: 873F 3'-ATTggCAAAGATAgCAggg-5', 873R 3'-AgCATCAATCggCTgTggA-5', 884F 3'-ggTCAGCCTgTAAAAACACCG-5' and 884R 3'-gTCCCATCCATTAACCTTCA-5'; and for T5S-T4: 2341F 3'-ACACCACATCATTCACTggAgA-5', 2341R 3'-gCCTTCTAACATCCTgCA-5', 3296F 3'-gCTACATATCCAACCTCCTAgA-5' and 3296R 3'-AgCAATAgTAACCGACTTCTCAg-5' were used. Strain-specific primers were tested in 20 µL standard PCR reactions. The reaction mixture contained 1X Taq DNA polymerase buffer (10X: 500 mM KCl, 100 mM Tris-HCl, pH 9.0 at 25°C, 1% Triton®X-100), 200 µM of each dNTP, 10 pmol of each primer, and 0.5 U of the enzyme Taq DNA polymerase (Inbio Highway, Argentina). Amplifications were performed in a thermal cycler multigene TM II (Labnet International Inc., USA). During the first cycle, DNA was denatured at 94 °C for 5 min. For the subsequent 30 cycles, the tubes were kept at 94 °C for 40 s, at 57 °C for 70 s and at 72 °C for 65 s, followed by an extension at 72 °C for 10 min. The annealing temperatures of 55; 57, and 60 °C were tested. The PCR products were visualized in 2% (w/v) agarose gel stained with GelRed® (Sigma-Aldrich, Germany). The electrophoretic run was performed in an electrophoretic vessel (Electrophoresis Subsystem 70 Labnet International) at 110 V for 30 min and bands were visualized using a UV transilluminator (Model MUV 21-312-220). **Results:** The primers were specific for each *B. altitudinis* strain, as no amplification products were obtained for negative controls. Amplification products were produced using an annealing temperature of 57 °C. Amplicons of approximately 500

¹ Instituto de Biotecnología Misiones "Dra. María Ebe Reca" (InBioMis)- CONICET. Facultad de Ciencias Exactas Químicas y Naturales/FCEQyN. Universidad Nacional de Misiones/UNaM, cortesejulietta@gmail.com

² Instituto de Biotecnología Misiones "Dra. María Ebe Reca" (InBioMis), lolicastrillo82@gmail.com

³ Instituto de Biotecnología Misiones "Dra. María Ebe Reca" (InBioMis)- CONICET. Facultad de Ciencias Exactas Químicas y Naturales/FCEQyN. Universidad Nacional de Misiones/UNaM, onettoandrea@gmail.com

⁴ Instituto de Biotecnología Misiones "Dra. María Ebe Reca" (InBioMis)- CONICET. Facultad de Ciencias Exactas Químicas y Naturales/FCEQyN. Universidad Nacional de Misiones/UNaM, pdr_dario@yahoo.com

⁵ Faculty of Exact Chemical and Natural Sciences (FCEQyN) - National University of Misiones (UNaM), melaczkeski@fceqyn.unam.edu.ar

and 300 bp were generated for *B. altitudinis* 19RS3 and T5S-T4, respectively. Conclusions: *B. altitudinis* 19RS3 and *B. altitudinis* T5S-T4 were successfully detected. The use of strain-specific primers are one of the cheapest and quickest methods to monitoring the colonization of bacterial strains in nursery and field experiments. The strain-specific primers will be applied in future traceability experiments.

PALAVRAS-CHAVE: Biofertilizer, Degenerated primers, Monitoring, Plant growth promoting bacteria, Traceability

¹ Instituto de Biotecnología Misiones "Dra. María Ebe Reca" (InBioMis)- CONICET. Facultad de Ciencias Exactas Químicas y Naturales/FCEQyN. Universidad Nacional de Misiones/UNaM, cortesejulietta@gmail.com
² Instituto de Biotecnología Misiones "Dra. María Ebe Reca" (InBioMis), lolicastrillo82@gmail.com
³ Instituto de Biotecnología Misiones "Dra. María Ebe Reca" (InBioMis)- CONICET. Facultad de Ciencias Exactas Químicas y Naturales/FCEQyN. Universidad Nacional de Misiones/UNaM, onettoandrea@gmail.com
⁴ Instituto de Biotecnología Misiones "Dra. María Ebe Reca" (InBioMis)- CONICET. Facultad de Ciencias Exactas Químicas y Naturales/FCEQyN. Universidad Nacional de Misiones/UNaM, pdr_dario@yahoo.com
⁵ Faculty of Exact Chemical and Natural Sciences (FCEQyN) - National University of Misiones (UNaM)., melaczeski@fceqyn.unam.edu.ar